

Amersham beads PCR Polymerase-Chain-Reaction



Procedure 1

using 0,2 ml-Tubes
PCR-machine with
lid-heating



Variant A

using
Amersham-Beads
with
Buffer/Polymerase



Batch 1A

2,5 µl Template
1 µl Primer forward
1 µl Primer reverse
20 µl dd H₂O



Batch 2A

= 2x Batch 1A
place tubes on ice until
PCR starts



Procedure 2

using 0,5 ml-tubes
PCR-machine without
lid-heating



Variant B

classisc procedure
all reagents are
added seperately



Batch 1B

2,5 µl Template
4 µl dNTP-Mix
13 µl dd H₂O
1 µl Taq-Polymerase
1 µl Primer forward
1 µl Primer reverse



Batch 2B

= 2x Batch 1B
place tubes on ice until
PCR starts




optional:

add running buffer
Kresolrot in 6M Glc



Program PCR-
machine as advised
or
call available program



95°C 60s Lead in
Cycle 1-35
XX°C 30s Annealing
72°C 90s Polymerize
95°C 30s Denature



XX°C 30s Annealing
72°C 3m Polymerize
95°C 5m Leadout
4°C store
XX= according to primer