

## extract genomic DNA



Prepare 3 ml H<sub>2</sub>O or 0,9% NaCl-solution for collecting mucosa-cells



Chew inner cheeks gently and flush liquid back into the falcon tube.




Mix the cell-suspension well for 30 sec on a vortex to get a homogen solution



Transfer 1ml for into Epi-tube for subsequent DNA - extraction



Sediment cells through short spinning for 30 s decant supernatant and add 30 µl dd H<sub>2</sub>O



Add Chelex(R)-beads to remove DNase-activating metal-ions




Add 100 µl Chelex(R) and resuspend pellet by vortexing



Cell-lysis is achieved by heat shock for 10 min at 99 °C Use waterbath or metal-block heater



Sediment Chelex beads by short spinning for 30 sec Never miss balancing the rotor by using exact counterweight



Remove supernatant carefully.

Don't carry Chelex beads!



Supernatant contains the DNA. To stock at deep freeze transfer into an Epi-tube or use directly.



E.g. as template in a PCR-reaction. Take 2,5µl DNA + 22,5 µl PCR-Mix or H<sub>2</sub>O, if using Amersham PCR-Beads